

# MRI MULTI-DOSING

*Medical Systems sets a new standard of excellence in MRI multi-dosing with "Easy Connect™", an all-in-one automatic dual check valve kit for dual head injectors...*

**Easy Connect™**  
An Integrated MRI Multi-dosing Kit for Dual Head Injectors

**TR-60 2B/C**  
Patient Line with 2 Check Valves

**TR-8000**  
Saline Line

*A simple line change is all you need for each patient...*

Spike for Contrast Vial

1. Quick, simple & safe
2. 900 psi and 8cc/sec
3. No loose connections, cross threading or leakages
4. Less prone to touch contamination
5. Less Inventory to carry
6. Spike for Contrast Vial
7. Automatic Dual Check Valve in each line

**Essential Requirement**  
Health Canada's preliminary report of February 15, 1996, on infection control, suggests that multiple patient dosing of expensive contrast media for achieving cost reductions is a safe option when transfer devices with appropriately designed and tested check valves are used.

Medical Systems' check valves and devices have been validated by an independent Canadian laboratory for preventing the backflow of potentially contaminated fluids in a multi-dosing setting of contrast media. Complete line of testing includes bacteriological, radiopharmacy, and structural testing to ensure the highest margin of safety on our check valves.

**MRI KIT for Liebel-Flarsheim**

**TR-801021 LF**

- 2 - 60ml Tyco Syringes
- 2 - Filling Spikes
- Y Tubing 60" with Check Valve

**High Flow Patient Line Extension sets with 2 check valves**

- TR-60 2B/C - 60"
- TR-72 2B/C - 72"
- TR-05 2B/C - 5"
- SR-40 2B/C - 40"

- High pressure 500 psi lines
- High flow 8cc/ sec
- Superior flexibility
- No blown check valves
- No blown tubing
- No leakages

**ALL PRODUCTS ARE LATEX FREE & NON-DEHP**

# Medical Systems



**We take your business To Heart...**

P.O. Box 48104, Mississauga, Ont L5A 4G8 Canada  
Tel: 1-888-344 4491 Fax: 1-866-553 5508  
Email: [medsystems@rogers.com](mailto:medsystems@rogers.com) [www.medicalsystems.ca](http://www.medicalsystems.ca)

Brochure-MRI Multidosing Rev. 1 05/20/2008

Please contact Stephen Romaniuk, National Sales Manager, for all your enquiries or questions - toll free at 1-866-636 2667