

KIT PACKAGING



Let Medical Systems help you to customize your kit requirements!!

We will customize your kits to any specifications you desire.....with components and/or finished medical devices to save you on inventory hassles and cost..... for General Surgery, O.R., Vein Clinics, Cardiac Cath Lab, Angiography Department and others.....



Equipment Covers



Stopcocks, Manifolds



Guidewires




I.V. Sets




Pressure Lines up to 1200 PSI




Waste Disposal System



Angio Drapes



Fluoroscope Drapes



Probe Covers



Gowns



Half Drapes



Surgical Gloves

Nerve Root Block Kit

- 1 Table Cover 44" x 90"
- 1 Banded Bag, 30" x 36"
- 2 Dome Bag, Round 40"
- 1 Syringe 5ml, L.L.
- 1 Syringe 10ml, L.L.
- 1 Control Syringe 10ml L.L.
- 1 16" Extension Set w/ Male/Female LL Connector & Roller Clamp



- 1 Needle 25G x 1.5"
- 3 Needle 18G x 1.5"

Arterial Kit

- 1 Table Cover 44" x 90"
- 1 Banded Bag, 30" x 36"
- 2 Dome Bag, Round 40"
- 1 Syringe 5ml, L.L.
- 4 Syringe 10ml, L.L.
- 1 Control Syringe 10ml L.L.
- 1 Telfa Pad 3" x 8"
- 1 Introducer Needle w/Base Plate
- 1 Waste Disposal System w/ Dual Check Valve, 2000ml waste Bag & 72" Fluid Administration Line
- Labels - 4 Blue (contrast), 5 Green (Hep N/S), 1 White, 1 Yellow & 3 Red



- 1 Needle 25G x 1.5"
- 1 Needle 18G x 1.5"
- 1 Scalpel #11

Varicocele Kit

- 1 Table Cover 44" x 90"
- 2 Banded Bag, 30" x 36"
- 2 Dome Bag, Round 40"
- 1 Fluoroscope Drape 35" x 43"
- 1 Syringe 20ml, L.L.
- 2 Syringe 10ml, L.L.
- 1 Control Syringe 10ml L.L.
- 1 Introducer Needle w/Base Plate
- 1 Probe Cover 6" w x 48"L w/5" cup, Gel & 2 Clips
- Labels - 2 Blue (contrast), 3 Green (Hep N/S), 3 White & 1 Red



- 1 Needle 25G x 1.5"
- 1 Needle 18G x 1.5"
- 1 Scalpel #11
- 1 Telfa Pad 3" x 8"

Vein (EVLT) Kit

- 1 Angio Drape 75" x 125"
- 2 Half Drape 44" x 55"
- 2 Gown Large
- 1 Surgical Glove # 6-1/2
- 1 Surgical Glove # 7
- 2 White Absorbent Towel
- 3 Steristrip 1/2" x 4"
- 1 Surgical Marker with Ruler
- 1 Probe Cover 6"w x 48"L w/5" cup & Gel
- 1 Needle 25G x 1.5"
- 3 Needle 22G x 1.5"
- 3 Syringe 20ml, L.L.
- 2 Syringe 10ml, L.L.
- 1 Blade #11 S.S.
- 2 Bowl 500cc
- 20 Gauze 4 x 4



Angiographic Kit

- 2 Banded Bag, 30" x 36"
- 1 Dome Bag, Round 40"
- 2 Syringe 20ml, L.L.
- 1 Syringe 20ml, Luer Slip
- 1 Syringe 10ml, L.L. (XYLOCAINE)
- 1 Syringe 10ml, Luer Slip (NITRO)
- 1 Syringe 3ml, L.L.
- 1 Syringe 1ml, Tuberculin Slip Tip
- 1 Braided 1200 PSI Extension Set, 30", w/ Male & Female Luer Lock
- 1 Manifold 3 Gang with (1) Waste Disposal System w/ high pressure automatic dual check valve, waste bag and Fluid Administration Line with infusion chamber, (2) Fluid Administration Line with vented spike, roller clamp and male luer lock, and (3) Pressure Line with male and female luer lock connectors
- Labels - 2 White (VERAPAMIL), 2 Red (NITRO), 1 Yellow, 1 Blue & 1 Red
- 1 Needle 16G x 1.5"
- 1 Needle 18G x 1.5"
- 1 Needle 20G x 1.5"
- 1 Needle 25G x 1.5"
- 1 Needle Holder
- 1 Scalpel #11
- 1 Control Syringe 10ml L.L.
- 1 Syringe Cannula
- 1 Closed Male Luer Lock Cap



Medical Systems

We take your business To Heart...

P.O. Box 48104, Mississauga, Ont L5A 4G8 Canada
 Tel: 1-888-344 4491 Fax: 1-866-553 5508
medsystems@rogers.com www.medicalsystems.ca



Made in Canada

ALL PRODUCTS ARE LATEX FREE & NON-DEHP

Please contact Stephen Romaniuk, National Sales Manager, for all your enquiries or questions - toll free at 1-866-636 2667