

DISPENSING PINS

Large Bore Rapid Flow Vented Spike w/ 3 micron antibacterial hydrophobic air filter



Capped female luer lock for syringe deployment

902200

Mini Vented Dispensing Pin with 3 micron hydrophobic filter for drug vial reconstitution



Capped female luer lock for syringe deployment

MS-PH-111

Micro Vented Dispensing Pin with 3 micron hydrophobic filter for drug vial reconstitution



Capped female luer lock for syringe deployment

MS-PH-112

Non-vented Dispensing Pin with Needleless two-way valve



MS-3500LDP

Male luer lock syringe access for aspiration or injection of medication from inverted bags

MS-2500LDP (same as above but vented spike)

One Way Fluid Transfer Valve, Vented Spike with Hinged cap, for withdrawal of diluent or injectable Medication into syringes



One-way valve eliminates potential for fluid to re-enter the IV bag

Hinged Closure Cap

MS-3005

Chemotherapy Reconstitution Spikes with 0.2 micron hydrophobic filter



901100

901100-NL

Needleless Luer connector for closed system multi-dose access

Taky / Chemo Spike with 0.2 micron large surface area hydrophobic filter



Larger Venting Surface for faster withdrawal and injection

903300

Medical Systems

All Products are
LATEX FREE & NON-DEHP



Made in Canada

We take your business To Heart...

P.O. Box 48104, Mississauga, Ont L5A 4G8 Canada

Tel: 1-888-344 4491 Fax: 1-866-553 5508

Email: medsystems@rogers.com www.medicalsystems.ca

Brochure-Dispensing Pins Rev. 0 05/12/2008

Please contact Stephen Romaniuk, National Sales Manager,
for all your enquiries or questions - toll free at 1-866-636 2667